

Integrin alpha V (Capture)

DG-461C

Package : 可询价

Product Name : Integrin alpha V (Capture)**Cat.No.:** DG-461C**Application** : ELISA, CMIA, CLIA, CBA**Reactivity** : Human**Host species** : Rabbit

| | |
|--------------------------|---|
| Gene ID | 3685 |
| Swiss Prot | P06756 |
| Reactivity | Human |
| Application | ELISA, CMIA, CLIA, CBA |
| Host species | Rabbit |
| Clonality | Monoclonal |
| Isotype | IgG |
| Purity | Affinity Purification |
| Concentration | 2 mg/ml |
| Storage Stability | Store at -20°C. Recommended to aliquot into single-use vials. Supplied in 1X PBS (pH 7.4). BSA and Azide Free. Stable for 12 months from date of receipt. |