

Growth hormone receptor (Capture)**DG-723C****Package : 可询价****Product Name :** Growth hormone receptor (Capture)**Cat.No.:** DG-723C**Application :** ELISA, CMIA, CLIA, CBA**Reactivity :** Human**Host species :** Rabbit**Gene ID** 2690**Swiss Prot** P10912**Reactivity** Human**Application** ELISA, CMIA, CLIA, CBA**Host species** Rabbit**Clonality** Monoclonal**Isotype** IgG**Purity** Affinity Purification**Concentration** 2mg/mL**Storage Stability** Store at -20°C. Recommended to aliquot into single-use vials. Supplied in 1X PBS (pH 7.4). BSA and Azide Free. Stable for 12 months from date of receipt.