

## Growth Hormone (Capture)

DG-7C

Package : 可询价

**Product Name :** Growth Hormone (Capture)

**Cat.No.:** DG-7C

**Application :** ELISA, CMIA, CLIA, CBA

**Reactivity :** Human

**Host species :** Rabbit

**Gene ID** 2688

**Swiss Prot** P01241

**Reactivity** Human

**Application** ELISA, CMIA, CLIA, CBA

**Host species** Rabbit

**Clonality** Monoclonal

**Isotype** IgG

**Purity** Affinity Purification

**Concentration** 2 mg/ml

**Storage Stability** Store at -20°C. Recommended to aliquot into single-use vials. Supplied in 1X PBS (pH 7.4). BSA and Azide Free. Stable for 12 months from date of receipt.