

Vitamin D Binding protein (Capture)

DG-334C

Package : 可询价

Product Name : Vitamin D Binding protein (Capture)**Cat.No.:** DG-334C**Application** : ELISA, CMIA, CLIA, CBA**Reactivity** : Human**Host species** : Rabbit

Gene ID	2638
Swiss Prot	P02774
Reactivity	Human
Application	ELISA, CMIA, CLIA, CBA
Host species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Purity	Affinity Purification
Concentration	2 mg/ml
Storage Stability	Store at -20°C. Recommended to aliquot into single-use vials. Supplied in 1X PBS (pH 7.4). BSA and Azide Free. Stable for 12 months from date of receipt.