

IL-2RA (Capture)**DG-1321C****Package : 可询价****Product Name** : IL-2RA (Capture)**Cat.No.:** DG-1321C**Application** : ELISA, CMIA, CLIA, CBA**Reactivity** : Human**Host species** : Rabbit

Gene ID	3559
Swiss Prot	P01589
Reactivity	Human
Application	ELISA, CMIA, CLIA, CBA
Host species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Purity	Affinity Purification
Concentration	2mg/mL
Storage Stability	Store at -20°C. Recommended to aliquot into single-use vials. Supplied in 1X PBS (pH 7.4). BSA and Azide Free. Stable for 12 months from date of receipt.