

**PCT (Detector)**

DG-279D

Package : 可询价

**Product Name** : PCT (Detector)**Cat.No.:** DG-279D**Application** : ELISA, CMIA, CLIA, CBA**Reactivity** : Human**Host species** : Rabbit**Gene ID** 796**Swiss Prot** P01258**Reactivity** Human**Application** ELISA, CMIA, CLIA, CBA**Host species** Rabbit**Clonality** Monoclonal**Isotype** IgG**Purity** Affinity Purification**Concentration** 2 mg/ml**Storage Stability** Store at -20°C. Recommended to aliquot into single-use vials. Supplied in 1X PBS (pH 7.4). BSA and Azide Free. Stable for 12 months from date of receipt.