

Phospho-Tau (Thr231) (Capture)**DG-1342C****Package : 可询价****Product Name** : Phospho-Tau (Thr231) (Capture)**Cat.No.:** DG-1342C**Application** : ELISA, CMIA, CLIA, CBA**Reactivity** : Human**Host species** : Rabbit**Gene ID** 4137**Swiss Prot** P10636**Reactivity** Human**Application** ELISA, CMIA, CLIA, CBA**Host species** Rabbit**Clonality** Monoclonal**Isotype** IgG**Purity** Affinity Purification**Concentration** 2 mg/ml**Storage Stability** Store at -20°C. Recommended to aliquot into single-use vials. Supplied in 1X PBS (pH 7.4). BSA and Azide Free. Stable for 12 months from date of receipt.