

ANGPTL3 (Capture)**DG-678C****Package : 可询价****Product Name :** ANGPTL3 (Capture)**Cat.No.:** DG-678C**Application :** ELISA, CMIA, CLIA, CBA**Reactivity :** Human**Host species :** Rabbit**Gene ID** 27329**Swiss Prot** Q9Y5C1**Reactivity** Human**Application** ELISA, CMIA, CLIA, CBA**Host species** Rabbit**Clonality** Monoclonal**Isotype** IgG**Purity** Affinity Purification**Concentration** 2 mg/ml**Storage Stability** Store at -20°C. Recommended to aliquot into single-use vials. Supplied in 1X PBS (pH 7.4). BSA and Azide Free. Stable for 12 months from date of receipt.