

ANGPTL3 (Detector)

DG-678D

Package : 可询价

Product Name : ANGPTL3 (Detector)**Cat.No.:** DG-678D**Application :** ELISA, CMIA, CLIA, CBA**Reactivity :** Human**Host species :** Rabbit**Gene ID** 27329**Swiss Prot** Q9Y5C1**Reactivity** Human**Application** ELISA, CMIA, CLIA, CBA**Host species** Rabbit**Clonality** Monoclonal**Isotype** IgG**Purity** Affinity Purification**Concentration** 2 mg/ml**Storage Stability** Store at -20°C. Recommended to aliquot into single-use vials. Supplied in 1X PBS (pH 7.4). BSA and Azide Free. Stable for 12 months from date of receipt.