

Recombinant Human Visfatin Protein (Standard)

DG-768P

Package : 可询价

Product Name : Recombinant Human Visfatin Protein (Standard)

Cat.No.: DG-768P

Synonyms : Pre-B-cell colony-enhancing factor 1, NAMPT, MGC117256, Nicotinamide phosphoribosyltransferase precursor, PBEF, PBEF1, VF, NAMPT, Visfatin, Nicotinamide phosphoribosyltransferase precursor, pre-B-cell colony enhancing factor 1,

Application : ELISA

Expression host : E.coli

Species : Human

Immunogen	Recombinant protein
Gene ID	10135
Swiss Prot	P43490
Synonyms	Pre-B-cell colony-enhancing factor 1, NAMPT, MGC117256, Nicotinamide phosphoribosyltransferase precursor, PBEF, PBEF1, VF, NAMPT, Visfatin, Nicotinamide phosphoribosyltransferase precursor, pre-B-cell colony enhancing factor 1,
Application	ELISA
Molecular Weight	57 kDa (511 aa), (SDS-PAGE under reducing conditions).
Species	Human
Purity	>90% as determined by SDS-PAGE.
Seq	1-491aa
Concentration	1 mg/ml (determined by Bradford assay)
Formulation	Freeze-dried powder
Detection Method	Sandwich ELISA
Storage Stability	Store at 2-8℃ for one month. Aliquot and store at -80℃ for 12 months. Avoid repeated freeze/thaw cycles.

□ 15% SDS PAGE(reduced) analysis of Recombinant Human Visfatin Protein (Standard)