

Recombinant Human C-X-C Motif Chemokine 7/CXCL7/NAP-2 Protein (Standard)

DG-222P

Package : 可询价

Product Name : Recombinant Human C-X-C Motif Chemokine 7/CXCL7/NAP-2 Protein (Standard)**Cat.No.**: DG-222P**Synonyms** : Platelet Basic Protein; PBP; C-X-C Motif Chemokine 7; Leukocyte-Derived Growth Factor; LDGF; Macrophage-Derived Growth Factor; MDGFSmall-Inducible Cytokine B7; PPBP; CTAP3; CXCL7; SCYB7; TGB1; THBGB1**Application** : ELISA**Expression host** : E.coli**Species** : Human

Immunogen	Recombinant protein
Gene ID	5473
Swiss Prot	P02775
Synonyms	Platelet Basic Protein; PBP; C-X-C Motif Chemokine 7; Leukocyte-Derived Growth Factor; LDGF; Macrophage-Derived Growth Factor; MDGFSmall-Inducible Cytokine B7; PPBP; CTAP3; CXCL7; SCYB7; TGB1; THBGB1
Application	ELISA
Molecular Weight	7.6 kDa
Species	Human
Purity	>90% as determined by SDS-PAGE.
Seq	Ala59-Asp128
Formulation	Freeze-dried powder
Detection Method	Sandwich ELISA
Storage Stability	Store at 2-8℃ for one month. Aliquot and store at -80℃ for 12 months. Avoid repeated freeze/thaw cycles.
	15% SDS PAGE(reduced) analysis of Recombinant Human C-X-C Motif Chemokine 7/CXCL7/NAP-2 Protein (Standard)