

**Alpha 2 Antiplasmin (Detector)****DG-670D****Package : 可询价****Product Name** : Alpha 2 Antiplasmin (Detector)**Cat.No.:** DG-670D**Application** : ELISA, CMIA, CLIA, CBA**Reactivity** : Human**Host species** : Rabbit

|                          |                                                                                                                                                           |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene ID</b>           | 5345                                                                                                                                                      |
| <b>Swiss Prot</b>        | P08697                                                                                                                                                    |
| <b>Reactivity</b>        | Human                                                                                                                                                     |
| <b>Application</b>       | ELISA, CMIA, CLIA, CBA                                                                                                                                    |
| <b>Host species</b>      | Rabbit                                                                                                                                                    |
| <b>Clonality</b>         | Monoclonal                                                                                                                                                |
| <b>Isotype</b>           | IgG                                                                                                                                                       |
| <b>Purity</b>            | Affinity Purification                                                                                                                                     |
| <b>Concentration</b>     | 2 mg/ml                                                                                                                                                   |
| <b>Storage Stability</b> | Store at -20°C. Recommended to aliquot into single-use vials. Supplied in 1X PBS (pH 7.4). BSA and Azide Free. Stable for 12 months from date of receipt. |