

CD276 Rabbit mAb

PLD-630

Package :

Product Name : CD276 Rabbit mAb**Cat.No.:** PLD-630**Application** : IHC-P**Reactivity** : Human**Reactivity** Human**Application** IHC-P**Recommended dilution** IHC-P: 1:200**Concentration** 2mg/mL**Storage Stability** Store at -20°C. Recommended to aliquot into single-use vials. Supplied in 1X PBS (pH 7.4). BSA and Azide Free. Stable for 12 months from date of receipt.