

MMP9 Rabbit mAb**PLD-115****Package :****Product Name :** MMP9 Rabbit mAb**Cat.No.:** PLD-115**Application :** IHC-P**Reactivity :** Human**Reactivity** Human**Application** IHC-P**Recommended dilution** IHC-P: 1:200**Concentration** 2mg/mL**Storage Stability** Store at -20°C. Recommended to aliquot into single-use vials. Supplied in 1X PBS (pH 7.4). BSA and Azide Free. Stable for 12 months from date of receipt.