

**Prostate Specific Antigen (Capture)****DG-527C****Package : 可询价****Product Name :** Prostate Specific Antigen (Capture)**Cat.No.:** DG-527C**Application :** ELISA, CMIA, CLIA, CBA**Reactivity :** Human**Host species :** Rabbit

|                          |   |
|--------------------------|---|
| <b>Gene ID</b>           | 354   |
| <b>Swiss Prot</b>        | P07288  |
| <b>Reactivity</b>        | Human   |
| <b>Application</b>       | ELISA, CMIA, CLIA, CBA  |
| <b>Host species</b>      | Rabbit  |
| <b>Clonality</b>         | Monoclonal  |
| <b>Isotype</b>           | IgG   |
| <b>Purity</b>            | Affinity Purification   |
| <b>Concentration</b>     | 2mg/mL  |
| <b>Storage Stability</b> | Store at -20°C. Recommended to aliquot into single-use vials. Supplied in 1X PBS (pH 7.4). BSA and Azide Free. Stable for 12 months from date of receipt. |