

CEACAM6 (Capture)

DG-1040C

Package : 可询价

Product Name : CEACAM6 (Capture)

Cat.No.: DG-1040C

Application : ELISA, CMIA, CLIA, CBA

Reactivity : Human

Host species : Rabbit

Gene ID 4680

Swiss Prot P40199

Reactivity Human

Application ELISA, CMIA, CLIA, CBA

Host species Rabbit

Clonality Monoclonal

Isotype IgG

Purity Affinity Purification

Concentration 2 mg/ml

Storage Stability Store at -20°C. Recommended to aliquot into single-use vials. Supplied in 1X PBS (pH 7.4). BSA and Azide Free. Stable for 12 months from date of receipt.