

S100 β (Detector)

DG-518D

Package : 可询价

Product Name : S100 β (Detector)
Cat.No.: DG-518D
Application : ELISA, CMIA, CLIA, CBA
Reactivity : Human
Host species : Rabbit

Gene ID	6285
Swiss Prot	P04271
Reactivity	Human
Application	ELISA, CMIA, CLIA, CBA
Host species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Purity	Affinity Purification
Concentration	2 mg/ml
Storage Stability	Store at -20°C. Recommended to aliquot into single-use vials. Supplied in 1X PBS (pH 7.4). BSA and Azide Free. Stable for 12 months from date of receipt.