

IL-1 β (Capture)

DG-160C

Package : 可询价

Product Name : IL-1 β (Capture)

Cat.No.: DG-160C

Application : ELISA, CMIA, CLIA, CBA

Reactivity : Human

Host species : Rabbit

Gene ID 3553

Swiss Prot P01584

Reactivity Human

Application ELISA, CMIA, CLIA, CBA

Host species Rabbit

Clonality Monoclonal

Isotype IgG

Purity Affinity Purification

Concentration 2mg/mL

Storage Stability Store at -20°C. Recommended to aliquot into single-use vials. Supplied in 1X PBS (pH 7.4). BSA and Azide Free. Stable for 12 months from date of receipt.

- ELISA法检测结果
- 化学发免疫法检测结果