

**Tissue Polypeptide Specific Antigen (Capture)**

DG-683C

Package : 可询价

**Product Name** : Tissue Polypeptide Specific Antigen (Capture)**Cat.No.:** DG-683C**Application** : ELISA, CMIA, CLIA, CBA**Reactivity** : Human**Host species** : Rabbit**Gene ID** 3875**Swiss Prot** P05783**Reactivity** Human**Application** ELISA, CMIA, CLIA, CBA**Host species** Rabbit**Clonality** Monoclonal**Isotype** IgG**Purity** Affinity Purification**Concentration** 2 mg/ml**Storage Stability** Store at -20°C. Recommended to aliquot into single-use vials. Supplied in 1X PBS (pH 7.4). BSA and Azide Free. Stable for 12 months from date of receipt.