

Alkaline Phosphatase (Detector)

DG-671D

Package : 可询价

Product Name : Alkaline Phosphatase (Detector)

Cat.No.: DG-671D

Application : ELISA, CMIA, CLIA, CBA

Reactivity : Human

Host species : Rabbit

Gene ID 248、249、250、251

Swiss Prot P05186、P05187、P09923、P10696

Reactivity Human

Application ELISA, CMIA, CLIA, CBA

Host species Rabbit

Clonality Monoclonal

Isotype IgG

Purity Affinity Purification

Concentration 2mg/mL

Storage Stability Store at -20°C. Recommended to aliquot into single-use vials. Supplied in 1X PBS (pH 7.4). BSA and Azide Free. Stable for 12 months from date of receipt.