

## Transferrin (Capture)

DG-1C

Package : 可询价

**Product Name** : Transferrin (Capture)**Cat.No.:** DG-1C**Application** : ELISA, CMIA, CLIA, CBA**Reactivity** : Human**Host species** : Rabbit

|                   |                                                                                                                                                           |
|-------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------|
| Gene ID           | 7018                                                                                                                                                      |
| Swiss Prot        | P02787                                                                                                                                                    |
| Reactivity        | Human                                                                                                                                                     |
| Application       | ELISA, CMIA, CLIA, CBA                                                                                                                                    |
| Host species      | Rabbit                                                                                                                                                    |
| Clonality         | Monoclonal                                                                                                                                                |
| Isotype           | IgG                                                                                                                                                       |
| Purity            | Affinity Purification                                                                                                                                     |
| Concentration     | 2 mg/ml                                                                                                                                                   |
| Storage Stability | Store at -20°C. Recommended to aliquot into single-use vials. Supplied in 1X PBS (pH 7.4). BSA and Azide Free. Stable for 12 months from date of receipt. |