

Hu IgD (Human IgD heavy chain)

PLD-156

Package : 可询价

Product Name : Hu IgD (Human IgD heavy chain)**Cat.No.:** PLD-156**Synonyms** : IGHD**Application** : IHC-P**Reactivity** : Human**Host species** : Rabbit

Background	Predicted to enable antigen binding activity and immunoglobulin receptor binding activity. Involved in positive regulation of interleukin-1 production. Located in blood microparticle and extracellular exosome. [provided by Alliance of Genome Resources, Jul 2025]
Immunogen	A synthetic peptide of human IgD
Gene ID	3495
Swiss Prot	P01880
Synonyms	IGHD
Reactivity	Human
Application	IHC-P
Recommended dilution	IHC-P: 1:2000-1:4000
Host species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Purity	Affinity Purification
Conjugation	Un-conjugated
Storage Stability	Store at -20°C. Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% BSA. Stable for 12 months from date of receipt.