

**BCL-2**

PLD-31

Package : 可询价

**Product Name** : BCL-2  
**Cat.No.:** PLD-31  
**Synonyms** : Bcl-2; PPP1R50  
**Application** : IHC-P  
**Reactivity** : Human  
**Host species** : Rabbit

<b>Background</b>	This gene encodes an integral outer mitochondrial membrane protein that blocks the apoptotic death of some cells such as lymphocytes. Constitutive expression of BCL2, such as in the case of translocation of BCL2 to Ig heavy chain locus, is thought to be the cause of follicular lymphoma. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Feb 2016]
<b>Gene ID</b>	596
<b>Swiss Prot</b>	P10415
<b>Synonyms</b>	Bcl-2; PPP1R50
<b>Reactivity</b>	Human
<b>Application</b>	IHC-P
<b>Recommended dilution</b>	1:500
<b>Host species</b>	Rabbit
<b>Clonality</b>	Monoclonal
<b>Isotype</b>	IgG
<b>Purity</b>	Affinity Purification
<b>Conjugation</b>	Un-conjugated
<b>Storage Stability</b>	Store at -20°C. Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% BSA. Stable for 12 months from date of receipt.