

FoxC1

PLD-461

Package : 可询价

Product Name : FoxC1**Cat.No.:** PLD-461**Synonyms** : ARA; IGDA; IHG1; ASGD3; FKHL7; IRID1; RIEG3; FREAC3; FREAC-3**Application** : IHC-P**Reactivity** : Human**Host species** : Rabbit**Background**

This gene belongs to the forkhead family of transcription factors which is characterized by a distinct DNA-binding forkhead domain. The specific function of this gene has not yet been determined; however, it has been shown to play a role in the regulation of embryonic and ocular development. Mutations in this gene cause various glaucoma phenotypes including primary congenital glaucoma, autosomal dominant iridogoniodysgenesis anomaly, and Axenfeld-Rieger anomaly. [provided by RefSeq, Jul 2008]

Immunogen

Recombinant protein of human FOXC1

Gene ID

2296

Swiss Prot

Q12948

Synonyms

ARA; IGDA; IHG1; ASGD3; FKHL7; IRID1; RIEG3; FREAC3; FREAC-3

Reactivity

Human

Application

IHC-P

Recommended dilution

IHC-P: 1:500

Host species

Rabbit

Clonality

Monoclonal

Isotype

IgG

Purity

Affinity Purification

Conjugation

Un-conjugated

Storage Stability

Store at -20°C. Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% BSA. Stable for 12 months from date of receipt.