

IGJ

PLD-677

Package : 可询价

Product Name : IGJ
Cat.No.: PLD-677
Synonyms : IGJ; JCH; IGCJ
Application : IHC-P
Reactivity : Human
Host species : Rabbit

Background	Serves to link two monomer units of either IgM or IgA
Immunogen	Recombinant protein of human IGJ
Gene ID	3512
Swiss Prot	P01591
Synonyms	IGJ; JCH; IGCJ
Reactivity	Human
Application	IHC-P
Recommended dilution	IHC-P: 1:200
Host species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Purity	Affinity Purification
Conjugation	Un-conjugated
Storage Stability	Store at -20°C. Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% BSA. Stable for 12 months from date of receipt.