

**CRX**

PLD-312

Package : 可询价

**Product Name** : CRX**Cat.No.:** PLD-312**Synonyms** : CRD; LCA7; OTX3; CORD2**Application** : IHC-P**Reactivity** : Human**Host species** : Rabbit**Background**

The protein encoded by this gene is a photoreceptor-specific transcription factor which plays a role in the differentiation of photoreceptor cells. This homeodomain protein is necessary for the maintenance of normal cone and rod function. Mutations in this gene are associated with photoreceptor degeneration, Leber congenital amaurosis type III and the autosomal dominant cone-rod dystrophy 2. Several alternatively spliced transcript variants of this gene have been described, but the full-length nature of some variants has not been determined. [provided by RefSeq, Jul 2008]

**Gene ID**

1406

**Swiss Prot**

O43186

**Synonyms**

CRD; LCA7; OTX3; CORD2

**Reactivity**

Human

**Application**

IHC-P

**Recommended dilution**

1:200-1:500

**Host species**

Rabbit

**Clonality**

Monoclonal

**Isotype**

IgG

**Purity**

Affinity Purification

**Conjugation**

Un-conjugated

**Storage Stability**

Store at -20°C. Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% BSA. Stable for 12 months from date of receipt.