

GFAP

PLD-845

Package : 可询价

Product Name : GFAP
Cat.No.: PLD-845
Synonyms : ALXDRD
Application : IHC-P
Reactivity : Human
Host species : Rabbit

Background

This gene encodes one of the major intermediate filament proteins of mature astrocytes. It is used as a marker to distinguish astrocytes from other glial cells during development. Mutations in this gene cause Alexander disease, a rare disorder of astrocytes in the central nervous system. Alternative splicing results in multiple transcript variants encoding distinct isoforms. [provided by RefSeq, Oct 2008]

Gene ID

2670

Swiss Prot

P14136

Synonyms

ALXDRD

Reactivity

Human

Application

IHC-P

Recommended dilution

1:500-1:2000

Host species

Rabbit

Clonality

Monoclonal

Isotype

IgG

Purity

Affinity Purification

Conjugation

Un-conjugated

Storage Stability

Store at -20°C. Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% BSA. Stable for 12 months from date of receipt.