

**FOXJ1****PLD-486****Package : 可询价****Product Name** : FOXJ1**Cat.No.:** PLD-486**Synonyms** : HFH4; HFH-4; CILD43; FKHL13**Application** : IHC-P**Reactivity** : Human**Host species** : Rabbit**Background**

This gene encodes a member of the forkhead family of transcription factors. Similar genes in zebrafish and mouse have been shown to regulate the transcription of genes that control the production of motile cilia. The mouse ortholog also functions in the determination of left-right asymmetry. Polymorphisms in this gene are associated with systemic lupus erythematosus and allergic rhinitis.[provided by RefSeq, Sep 2009]

**Immunogen**

Recombinant protein of human FOXJ

**Gene ID**

2302

**Swiss Prot**

Q92949

**Synonyms**

HFH4; HFH-4; CILD43; FKHL13

**Reactivity**

Human

**Application**

IHC-P

**Recommended dilution**

IHC-P: 1:1000-1:2000

**Host species**

Rabbit

**Clonality**

Monoclonal

**Isotype**

IgG

**Purity**

Affinity Purification

**Conjugation**

Un-conjugated

**Storage Stability**

Store at -20°C. Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% BSA. Stable for 12 months from date of receipt.