

Recombinant Tissue Polypeptide Specific Antigen Protein (Standard)

DG-683P

Package : 可询价

Product Name : Recombinant Tissue Polypeptide Specific Antigen Protein (Standard)**Cat.No.:** DG-683P**Synonyms** : KRT18; CYK18; K18; Keratin 18; Cell proliferation-inducing gene 46 protein; Keratin, type I cytoskel**Application** : ELISA**Tag** : N-terminal His-Tag and Sumo**Expression host** : E.coli**Species** : Human

Immunogen	Recombinant protein
Gene ID	3875
Swiss Prot	P05783
Synonyms	KRT18; CYK18; K18; Keratin 18; Cell proliferation-inducing gene 46 protein; Keratin, type I cytoskel
Application	ELISA
Molecular Weight	36kDa
Species	Human
Purity	>90% as determined by SDS-PAGE.
Seq	Ala239-Asp397
Tag	N-terminal His-Tag and Sumo
Concentration	0.5mg/ml (determined by Absorbance at 280nm)
Formulation	Freeze-dried powder
Detection Method	Sandwich ELISA
Storage Stability	Aliquot and store at -20°C to -80°C for up to 6 months. Phosphate buffered saline (pH7.4) containing 0.01% sarcosyl, 5%Trehalose; buffer containing 50% Glycerol is recommen.