

Polymeric immunoglobulin receptor/PIGR (Detector)**DG-1230D****Package : 可询价****Product Name :** Polymeric immunoglobulin receptor/PIGR (Detector)**Cat.No.:** DG-1230D**Application :** ELISA, CMIA, CLIA, CBA**Reactivity :** Human**Host species :** Rabbit

Gene ID	5284
Swiss Prot	P01833
Reactivity	Human
Application	ELISA, CMIA, CLIA, CBA
Host species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Purity	Affinity Purification
Concentration	2 mg/ml
Storage Stability	Store at -20°C. Recommended to aliquot into single-use vials. Supplied in 1X PBS (pH 7.4). BSA and Azide Free. Stable for 12 months from date of receipt.