

## Polymeric immunoglobulin receptor/PIGR (Detector)

DG-1230D

Package : 可询价

**Product Name** : Polymeric immunoglobulin receptor/PIGR (Detector)**Cat.No.:** DG-1230D**Application** : ELISA, CMIA, CLIA, CBA**Reactivity** : Human**Host species** : Rabbit

<b>Gene ID</b>	5284
<b>Swiss Prot</b>	P01833
<b>Reactivity</b>	Human
<b>Application</b>	ELISA, CMIA, CLIA, CBA
<b>Host species</b>	Rabbit
<b>Clonality</b>	Monoclonal
<b>Isotype</b>	IgG
<b>Purity</b>	Affinity Purification
<b>Concentration</b>	2 mg/ml
<b>Storage Stability</b>	Store at -20°C. Recommended to aliquot into single-use vials. Supplied in 1X PBS (pH 7.4). BSA and Azide Free. Stable for 12 months from date of receipt.