

CXCL7/PBP (Detector)

DG-222D

Package : 可询价

Product Name : CXCL7/PBP (Detector)

Cat.No.: DG-222D

Application : ELISA, CMIA, CLIA, CBA

Reactivity : Human

Host species : Rabbit

Gene ID 5473

Swiss Prot P02775

Reactivity Human

Application ELISA, CMIA, CLIA, CBA

Host species Rabbit

Clonality Monoclonal

Isotype IgG

Purity Affinity Purification

Concentration 2 mg/ml

Storage Stability Store at -20°C. Recommended to aliquot into single-use vials. Supplied in 1X PBS (pH 7.4). BSA and Azide Free. Stable for 12 months from date of receipt.