

**Recombinant Neuron-Specific Enolase Protein (Standard)**

DG-513P

Package : 可询价

**Product Name** : Recombinant Neuron-Specific Enolase Protein (Standard)**Cat.No.:** DG-513P**Synonyms** : Enolase 2 (gamma, neuronal)**Application** : ELISA**Expression host** : E.coli**Species** : Human

<b>Immunogen</b>	Recombinant protein
<b>Gene ID</b>	2026
<b>Swiss Prot</b>	P09104
<b>Synonyms</b>	Enolase 2 (gamma, neuronal)
<b>Application</b>	ELISA
<b>Molecular Weight</b>	47 kDa ( 434 aa)
<b>Species</b>	Human
<b>Purity</b>	>85% as determined by SDS-PAGE.
<b>Seq</b>	1-434aa
<b>Concentration</b>	1 mg/ml (determined by Bradford assay)
<b>Formulation</b>	Freeze-dried powder
<b>Detection Method</b>	Sandwich ELISA
<b>Storage Stability</b>	Store at 2-8°C for one month. Aliquot and store at -80°C for 12 months. Avoid repeated freeze/thaw cycles.