

**CD22****PLD-491****Package : 可询价**

**Product Name :** CD22  
**Cat.No.:** PLD-491  
**Synonyms :** SIGLEC2; SIGLEC-2  
**Application :** IHC-P  
**Reactivity :** Human  
**Host species :** Rabbit

|                             |                                                                                                                                                               |
|-----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Background</b>           | Mediates B-cell B-cell interactions                                                                                                                           |
| <b>Gene ID</b>              | 933                                                                                                                                                           |
| <b>Swiss Prot</b>           | P20273                                                                                                                                                        |
| <b>Synonyms</b>             | SIGLEC2; SIGLEC-2                                                                                                                                             |
| <b>Reactivity</b>           | Human                                                                                                                                                         |
| <b>Application</b>          | IHC-P                                                                                                                                                         |
| <b>Recommended dilution</b> | 1:2000-1:10000                                                                                                                                                |
| <b>Host species</b>         | Rabbit                                                                                                                                                        |
| <b>Clonality</b>            | Monoclonal                                                                                                                                                    |
| <b>Isotype</b>              | IgG                                                                                                                                                           |
| <b>Purity</b>               | Affinity Purification                                                                                                                                         |
| <b>Conjugation</b>          | Un-conjugated                                                                                                                                                 |
| <b>Storage Stability</b>    | Store at -20°C. Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% BSA. Stable for 12 months from date of receipt. |