

CD19

PLD-619

Package : 可询价

Product Name : CD19
Cat.No.: PLD-619
Synonyms : B4; CVID3
Application : IHC-P
Reactivity : Human
Host species : Rabbit

Background

Lymphocytes proliferate and differentiate in response to various concentrations of different antigens. The ability of the B cell to respond in a specific, yet sensitive manner to the various antigens is achieved with the use of low-affinity antigen receptors. This gene encodes a cell surface molecule which assembles with the antigen receptor of B lymphocytes in order to decrease the threshold for antigen receptor-dependent stimulation. [provided by RefSeq, Jul 2008]

Gene ID

930

Swiss Prot

P15391

Synonyms

B4; CVID3

Reactivity

Human

Application

IHC-P

Recommended dilution

1:2000-1:10000

Host species

Rabbit

Clonality

Monoclonal

Isotype

IgG

Purity

Affinity Purification

Conjugation

Un-conjugated

Storage Stability

Store at -20°C. Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% BSA. Stable for 12 months from date of receipt.