

IL-3 (Detector)**DG-315D****Package : 可询价****Product Name :** IL-3 (Detector)**Cat.No.:** DG-315D**Application :** ELISA, CMIA, CLIA, CBA**Reactivity :** Human**Host species :** Rabbit**Gene ID** 3562**Swiss Prot** P08700**Reactivity** Human**Application** ELISA, CMIA, CLIA, CBA**Host species** Rabbit**Clonality** Monoclonal**Isotype** IgG**Purity** Affinity Purification**Concentration** 2 mg/ml**Storage Stability** Store at -20°C. Recommended to aliquot into single-use vials. Supplied in 1X PBS (pH 7.4). BSA and Azide Free. Stable for 12 months from date of receipt.