

Ribonuclease 3/ECP (Detector)

DG-1245D

Package : 可询价

Product Name : Ribonuclease 3/ECP (Detector)

Cat.No.: DG-1245D

Application : ELISA, CMIA, CLIA, CBA

Reactivity : Human

Host species : Rabbit

Gene ID 6037

Swiss Prot P12724

Reactivity Human

Application ELISA, CMIA, CLIA, CBA

Host species Rabbit

Clonality Monoclonal

Isotype IgG

Purity Affinity Purification

Concentration 2 mg/ml

Storage Stability Store at -20°C. Recommended to aliquot into single-use vials. Supplied in 1X PBS (pH 7.4). BSA and Azide Free. Stable for 12 months from date of receipt.