

Amphiregulin (Capture)**DG-328C****Package : 可询价****Product Name :** Amphiregulin (Capture)**Cat.No.:** DG-328C**Application :** ELISA, CMIA, CLIA, CBA**Reactivity :** Human**Host species :** Rabbit

Gene ID	374
Swiss Prot	P15514
Reactivity	Human
Application	ELISA, CMIA, CLIA, CBA
Host species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Purity	Affinity Purification
Concentration	2 mg/ml
Storage Stability	Store at -20°C. Recommended to aliquot into single-use vials. Supplied in 1X PBS (pH 7.4). BSA and Azide Free. Stable for 12 months from date of receipt.