

**Glypican 1/ GPC1 (Capture)****DG-1327C****Package : 可询价****Product Name :** Glypican 1/ GPC1 (Capture)**Cat.No.:** DG-1327C**Application :** ELISA, CMIA, CLIA, CBA**Reactivity :** Human**Host species :** Rabbit

<b>Gene ID</b>	2817
<b>Swiss Prot</b>	P35052
<b>Reactivity</b>	Human
<b>Application</b>	ELISA, CMIA, CLIA, CBA
<b>Host species</b>	Rabbit
<b>Clonality</b>	Monoclonal
<b>Isotype</b>	IgG
<b>Purity</b>	Affinity Purification
<b>Concentration</b>	2 mg/ml
<b>Storage Stability</b>	Store at -20°C. Recommended to aliquot into single-use vials. Supplied in 1X PBS (pH 7.4). BSA and Azide Free. Stable for 12 months from date of receipt.