

Annexin-4/ANXA4 (Detector)

DG-960D

Package : 可询价

Product Name : Annexin-4/ANXA4 (Detector)

Cat.No.: DG-960D

Application : ELISA, CMIA, CLIA, CBA

Reactivity : Human

Host species : Rabbit

Gene ID	307
Swiss Prot	P09525
Reactivity	Human
Application	ELISA, CMIA, CLIA, CBA
Host species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Purity	Affinity Purification
Concentration	2 mg/ml
Storage Stability	Store at -20°C. Recommended to aliquot into single-use vials. Supplied in 1X PBS (pH 7.4). BSA and Azide Free. Stable for 12 months from date of receipt.