

CTLA-4 (Detector)

DG-762D

Package : 可询价

Product Name : CTLA-4 (Detector)**Cat.No.:** DG-762D**Application** : ELISA, CMIA, CLIA, CBA**Reactivity** : Human**Host species** : Rabbit**Gene ID** 1493**Swiss Prot** P16410**Reactivity** Human**Application** ELISA, CMIA, CLIA, CBA**Host species** Rabbit**Clonality** Monoclonal**Isotype** IgG**Purity** Affinity Purification**Concentration** 2 mg/ml**Storage Stability** Store at -20°C. Recommended to aliquot into single-use vials. Supplied in 1X PBS (pH 7.4). BSA and Azide Free. Stable for 12 months from date of receipt.