

Glutathione Peroxidase 3 (Detector)

DG-1373D

Package : 可询价

Product Name : Glutathione Peroxidase 3 (Detector)**Cat.No.:** DG-1373D**Application** : ELISA, CMIA, CLIA, CBA**Reactivity** : Human**Host species** : Rabbit**Gene ID** 2878**Swiss Prot** P22352**Reactivity** Human**Application** ELISA, CMIA, CLIA, CBA**Host species** Rabbit**Clonality** Monoclonal**Isotype** IgG**Purity** Affinity Purification**Concentration** 2 mg/ml**Storage Stability** Store at -20°C. Recommended to aliquot into single-use vials. Supplied in 1X PBS (pH 7.4). BSA and Azide Free. Stable for 12 months from date of receipt.