

14-3-3 sigma/SFN (Detector)**DG-588D****Package : 可询价****Product Name :** 14-3-3 sigma/SFN (Detector)**Cat.No.:** DG-588D**Application :** ELISA, CMIA, CLIA, CBA**Reactivity :** Human**Host species :** Rabbit

Gene ID	2810
Swiss Prot	P31947
Reactivity	Human
Application	ELISA, CMIA, CLIA, CBA
Host species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Purity	Affinity Purification
Concentration	2 mg/ml
Storage Stability	Store at -20°C. Recommended to aliquot into single-use vials. Supplied in 1X PBS (pH 7.4). BSA and Azide Free. Stable for 12 months from date of receipt.