

IGFBP4 (Detector)**DG-389D****Package : 可询价****Product Name :** IGFBP4 (Detector)**Cat.No.:** DG-389D**Application :** ELISA, CMIA, CLIA, CBA**Reactivity :** Human**Host species :** Rabbit

Gene ID	3487
Swiss Prot	P22692
Reactivity	Human
Application	ELISA, CMIA, CLIA, CBA
Host species	Rabbit
Clonality	Monoclonal Antibody
Isotype	IgG
Purity	Affinity Purification
Concentration	2 mg/ml
Storage Stability	Store at -20°C. Recommended to aliquot into single-use vials. Supplied in 1X PBS (pH 7.4). BSA and Azide Free. Stable for 12 months from date of receipt.