

IL-1 alpha (Detector)**DG-530D****Package : 可询价**

Product Name : IL-1 alpha (Detector)
Cat.No.: DG-530D
Application : ELISA, CMIA, CLIA, CBA
Reactivity : Human
Host species : Rabbit

Gene ID	3552
Swiss Prot	P01583
Reactivity	Human
Application	ELISA, CMIA, CLIA, CBA
Host species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Purity	Affinity Purification
Concentration	2 mg/ml
Storage Stability	Store at -20°C. Recommended to aliquot into single-use vials. Supplied in 1X PBS (pH 7.4). BSA and Azide Free. Stable for 12 months from date of receipt.