

**IL-1 $\beta$  (Capture)****DG-160C****Package : 可询价****Product Name :** IL-1 $\beta$  (Capture)**Cat.No.:** DG-160C**Application :** ELISA, CMIA, CLIA, CBA**Reactivity :** Human**Host species :** Rabbit

<b>Gene ID</b>	3553
<b>Swiss Prot</b>	P01584
<b>Reactivity</b>	Human
<b>Application</b>	ELISA, CMIA, CLIA, CBA
<b>Host species</b>	Rabbit
<b>Clonality</b>	Monoclonal
<b>Isotype</b>	IgG
<b>Purity</b>	Affinity Purification
<b>Concentration</b>	2 mg/ml
<b>Storage Stability</b>	Store at -20°C. Recommended to aliquot into single-use vials. Supplied in 1X PBS (pH 7.4). BSA and Azide Free. Stable for 12 months from date of receipt.
	<input type="checkbox"/> ELISA法检测结果
	<input type="checkbox"/> 化学发免疫法检测结果